




Enzymatic fingerprinting of Isomalto/Malto-polysaccharides

C.E. Klostermann*, P.H. van der Zaal*, H.A. Schols, J.H. Bitter & P.L. Buwalda

*these authors contributed equally to this work


Introduction

Starch consists of amylose & amylopectin, α -linked glucose polymers



Starch can be modified in different ways to increase its techno-functional properties. The enzymatic modification of starch with 4,6- α -glucanotransferase- Δ N (GTFB) results in Isomalto/Malto-polysaccharides (IMMPs). IMMPs are molecules built of α -1,4 linked glucose chains, connected to α -1,6 linked glucose chains (Leemhuis et al., 2014).

Graphical representation of enzymatic modification of starch by 4,6- α -glucanotransferase- Δ N, producing Isomalto/Malto-polysaccharides



Objective

"To develop an enzymatic fingerprinting method, capable to characterise linear Isomalto/Malto-polysaccharides"


Approach

In order to characterise linear IMMPs three enzymes were used:

- β -amylase for hydrolysis of all α -1,4 linkages from the non-reducing end
- Dextranase for the endo-hydrolysis of all α -1,6 linkages
- Isopullulanase to split hybrid α -1,4/ α -1,6 linked molecules


When using these enzymes in specific orders structural information can be obtained, as shown in the figure below.

Characterisation of linear IMMPs with β -amylase, dextranase and Isopullulanase. Enzymes with a "+" are incubated simultaneously, enzymes with a ">" are incubated after each other.




Results

In order to determine the presence of linear α -1,4 or linear α -1,6 linked chains in the IMMP mixture. HPAEC was used to compare α -1,4 and α -1,6 linked oligomer standards to an IMMP mixture containing 41 % α -1,6 linkages.




The presence of completely linear α -1,6 linked oligomers indicates that GTFB is able to use glucose as an acceptor.

In order to determine the length of the α -1,4 linked chains within the IMMPs, a linear IMMP (70 % α -1,6) was incubated with isopullulanase or β -amylase and isopullulanase simultaneously.



GTFB uses α -1,4 linked material \geq DP 6 as donor substrate and it uses α -1,4 linked material $<$ DP 6 as acceptor substrate.

In order to determine the length of the α -1,6 linked chains within IMMPs, a linear IMMP (70 % α -1,6) was incubated with dextranase or β -amylase and isopullulanase simultaneously. HPSEC was used to analyse these samples in the high Mw range.



GTFB is able to produce IMMPs with α -1,6 linked glucose chains up to 15 kDa.

Conclusions

The enzymatic fingerprinting method is able to characterise IMMPs.

GTFB produces:

- Fully linear α -1,6 linked chains
- IMMPs with α -1,4 linked glucose segments up to DP 5
- IMMPs with α -1,6 linked glucose segments up to 15 kDa

Acknowledgement

This project is jointly funded by TKI and AVEBE. We thank AVEBE for supplying starches and materials, the Rijksuniversiteit Groningen for supplying the GTFB- Δ N enzyme & analysis equipment, and the Carbohydrate Competence Centre for its support.